



Kombinationsmöglichkeiten Kryostate und Thermostate

	T100	T150	T200	T200K	TC150	TC250	TC250K
K10 Ca. 3 Liter	K10T100 -10...+100°C A-10857		K10T200 -10...+100°C A-10858				
K21 Ca. 8 Liter	K21T100 -20...+100°C A-10859		K21T200 A-10860	K21T200K -20...+200°C A-10861	K21TC150 -20...+100°C A-10862	K21TC250 -20...+100°C A-10863	K21TC250K -20...+200°C A-10864
K22 Ca. 10 Liter	K22T100 -20...+100°C A-10865		K22T200 -20...+100°C A-10866	K22T200K -20...+200°C A-10867	K22TC150 -20...+100°C A-10868	K22TC250 -20...+100°C A-10869	K22TC250K -20...+200°C A-10870
K31 Ca. 11 Liter	K31T100 -30...+100°C A-10871		K31T200 -30...+100°C A-10872	K31T200K -30...+200°C A-10873	K31TC150 -30...+100°C A-10874	K31TC250 -30...+100°C A-10875	K32TC250K -30...+200°C A-10876
K32 Ca. 15 Liter	K32T100 -30...+100°C A-10877		K32T200 -30...+100°C A-10878	K32T200K -30...+200°C A-10879	K32TC150 -30...+100°C A-10880	K32TC250 -30...+100°C A-10881	K32TC250K -30...+200°C A-10882
KT Ca. 14 Liter	KT100 0...+100°C A-10883		KT200 0...+100°C A-10884	KT200K 0...+200°C A-10885	KTTC150 0...+100°C A-10886	KTTC250 0...+100°C A-10887	KTTC250K 0...+200°C A-10888
KT10 Ca. 14 Liter	KT10T100 -10...+100°C A-10889		KT10T200 -10...+100°C A-10890	KT10T200K -10...+200°C A-10891	KT10TC150 -10...+100°C A-10892	KT10TC250 -10...+100°C A-10893	KT10TC250K -10...+200°C A-10894
KT30 Ca. 14 Liter		KT30T150 -30...+100°C A-10895	KT30T200 -30...+100°C A-10896	KT30T200K -30...+200°C A-10897	KT30TC150 -30...+100°C A-10898	KT30TC250 -30...+100°C A-10899	KT30TC250K -30...+200°C A-10900
KT40 Ca. 14 Liter		KT40T150 -40...+100°C A-10901	KT40T200 -40...+100°C A-10902	KT40T200K -40...+200°C A-10903	KT40TC150 -40...+100°C A-10904	KT40TC250 -40...+100°C A-10905	KT40TC250K -40...+200°C A-10906
KB10 Ca. 22 Liter	KB10T100 -10...+100°C A-10907		KB10T200 -10...+100°C A-10908	KB10T200K -10...+200°C A-10909	KB10TC150 -10...+100°C A-10910	KB10TC250 -10...+100°C A-10911	KB10TC250K -10...+200°C A-10912
KB30 Ca. 22 Liter		KB30T150 -30...+100°C A-10913	KB30T200 -30...+100°C A-10914	KB30T200K -30...+200°C A-10915	KB30TC150 -30...+100°C A-10916	KB30TC250 -30...+100°C A-10917	KB30TC250K -30...+200°C A-10918
KB40 Ca. 22 Liter		KB40T150 -40...+100°C A-10919	KB40T200 -40...+100°C A-10920	KB40T200K -40...+200°C A-10921	KB40TC150 -40...+100°C A-10922	KB40TC250 -40...+100°C A-10923	KB40TC250K -40...+200°C A-10924
KB55 Ca. 22 Liter		KB55T150 -55...+100°C A-10925	KB55T200 -55...+100°C A-10926	KB55T200K -55...+200°C A-10927	KB55TC150 -55...+100°C A-10928	KB55TC250 -55...+100°C A-10929	KB55TC250K -55...+200°C A-10930
KB75 Ca. 12 Liter		KB75T150 -75...+100°C A-10931	KB75T200 -75...+100°C A-10932	KB75T200K -75...+200°C A-10933	KB75TC150 -75...+100°C A-10934	KB75TC250 -75...+100°C A-10935	KB75TC250K -75...+200°C A-10936

Selbstverständlich sind von den genannten Größen abweichende Abmessung auf Anfrage realisierbar.